



Specification

8.55057.0000 Universal NovaTag™ resin Novabiochem®

Specification

Color (visual)	white to yellow to beige	
Appearance of substance (visual)	beads	
Loading (Photometric determination of the Fmoc-chromophore liberated upon treatment with DBU/DMF)	0.30 - 0.60	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g

The polymer matrix is copoly (styren- 1% DVB) 100 - 200 mesh

Dr. Jörg Bauer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.