



Specification

8.55070.0000 HMBA-PEGA resin Novabiochem®

Specification

Color (visual)	white to yellow to beige	
Appearance of substance (visual)	beads	
Loading (determined from the substitution of the Fmoc-Leu loaded resin)	0.20 - 0.50	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g
dry resin %	lot specific result	%

The resin is sold swollen in ethanol.

Dr. Jörg Bauer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.