



Specification

1.00103.0000 Amidosulfuric acid for analysis EMSURE®

	Spec. Values	
Assay (acidimetric)	≥ 99.0	%
Chloride (Cl)	≤ 100	ppm
Nitrate (NO ₃)	≤ 20	ppm
Sulfate (SO ₄)	≤ 500	ppm
Heavy metals (as Pb)	≤ 5	ppm
Fe (Iron)	≤ 5	ppm
Residue on ignition	≤ 200	ppm

Kristin Haufschildt

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.