



Specification

1.01024.0000 1-Propanol for liquid chromatography LiChrosolv®

Specification		
Purity (GC)	≥ 99.8	%
Identity (IR)	conforms	
Evaporation residue	≤ 2.0	mg/l
Water	≤ 0.02	%
Acidity	≤ 0.0002	meq/g
Alkalinity	≤ 0.0002	meq/g
Transmission (at 230 nm)	≥ 70	%
Transmission (at 240 nm)	≥ 80	%
Transmission (from 270 nm)	≥ 98	%
Filtered by 0.2 µm filter		

Dr. Ute Volkwein

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.