



# Specification

1.01115.0000 Ammonium acetate EMPLURA®

Spec. Values		
Assay (alkalimetric)	≥ 96.0	%
Identity	passes test	
pH-value (5 %; water)	6.0 - 7.5	
Chloride (Cl)	≤ 0.002	%
Sulfate (SO <sub>4</sub> )	≤ 0.01	%
Heavy metals (as Pb)	≤ 0.0005	%
Fe (Iron)	≤ 0.001	%
Residue on ignition (as sulfate)	≤ 0.02	%
Water	≤ 2.5	%

Claudia Wiegand

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.