



Specification

1.03819.0000 Iron for analysis reduced, particle size 10 µm EMSURE®

| Specification | | |
|---------------------------------------|----------|---|
| Assay (cerimetric) | ≥ 99.5 | % |
| Substances insoluble in sulfuric acid | ≤ 0.1 | % |
| Chloride (Cl) | ≤ 0.002 | % |
| Sulphide (S) | ≤ 0.002 | % |
| As (Arsenic) | ≤ 0.0002 | % |
| Cu (Copper) | ≤ 0.002 | % |
| Mn (Manganese) | ≤ 0.002 | % |
| Pb (Lead) | ≤ 0.002 | % |
| Zn (Zinc) | ≤ 0.002 | % |

Dr. Dimitrij Ryvlin

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.