



Specification

1.07687.0000 Sucrose (saccharose) for biochem. Reag. Ph Eur

Specification		
Identity (wet chemistry)	passes test	
Appearance of solution (50 %; water)	passes test	
Color number	≤ 45	
Spec. rotation (α 20/D; 26 %; water)	66.3 - 67.0	°
Electrical conductivity	≤ 35	μS/cm
Sulfite (as SO ₂)	≤ 10	ppm
Dextrins	passes test	
Reducing sugars	passes test	
Residual solvents (ICH Q3C)	excluded by production process	
Loss on Drying (105 °C)	≤ 0.1	%
Endotoxins	< 250	I.U./g

Julian Hofmann

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.