



Specification

1.07960.0000 Entellan™ rapid mounting medium for microscopy

Specification

Refractive index (n 20/D)	1.492 - 1.500	
Density (d 20 °C/ 4 °C)	0.925 - 0.935	
Fluorescence (as quinine at 365 nm)	≤ 100	ppb
Viscosity (20 °C)	60 - 100	mPa•s
Suitability for microscopy	passes test	

Dr. Timm Schlegelmilch

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.