



# Specification

---

1.09543.0000 pH-indicator strips, pH 6.5 - 10.0 non-bleeding pH 6.5 - 6.8 - 7.1 - 7.4 - 7.7 - 7.9 - 8.1 - 8.3 - 8.5 - 8.7 - 9.0 - 9.5 - 10.0 MQuant®

---

---

## Specification

---

Application test	conforms
------------------	----------

Sharon Brilmayer

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.