



Specification

1.10285.0000 LB broth (MILLER) for microbiology

Specification	
Appearance (clearness)	clear
Appearance (color)	yellowish to yellowish brown
pH-value (25 °C)	6.8 - 7.2

Specification	
Inoculum on reference medium (Escherichia coli ATCC 25922 (WDCM 00013))	≤ 100
Inoculum on reference medium (Escherichia coli ATCC 11775 (WDCM 00090))	≤ 100
Inoculum on reference medium (Escherichia coli DSM 6897 (DH5alpha))	≤ 100
Inoculum on reference medium (Escherichia coli ATCC 47092 (BL21))	≤ 100
Growth (Escherichia coli ATCC 25922 (WDCM 00013))	good to very good
Growth (Escherichia coli ATCC 11775 (WDCM 00090))	good to very good
Growth (Escherichia coli DSM 6897 (DH5alpha))	good to very good
Growth (Escherichia coli ATCC 47092 (BL21))	good to very good

Incubation: 24 ± 2 hours at 37 ± 1 °C, aerobic

Dr. Lukas Mechler

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.