



Specification

1.16738.0000 1-Methoxy-2-propanol EMPLURA®

Specification		
Purity (GC)	≥ 99.5	%
2-Methoxypropan-1-ol (GC)	≤ 0.29	%
Identity (IR)	conforms	
Free acid (as C ₂ H ₅ COOH)	≤ 0.002	%
Al (Aluminium)	≤ 20	ppb
Ca (Calcium)	≤ 50	ppb
Cu (Copper)	≤ 20	ppb
Fe (Iron)	≤ 100	ppb
K (Potassium)	≤ 50	ppb
Na (Sodium)	≤ 1000	ppb
Ni (Nickel)	≤ 20	ppb
Zn (Zinc)	≤ 20	ppb
Water	≤ 0.1	%

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.