



# Specification

---

1.16740.0000 Petroleum benzine boiling range 40-60°C hypergrade for organic trace analysis SupraSolv®

---

	Specification	
Boiling range (40-60°C)	≥ 90	% (v/v)
Evaporation residue	≤ 3.0	mg/l
Water	≤ 0.005	%
Color	≤ 10	Hazen
GC/ECD (retention range 1,2,4-tri-chlorobenzene to decachlorobiphenyle, individual signals (lindane standard))	≤ 2	pg/ml
GC/FID (retention range n-undecane to n-tetracontane, individual signals (n-tetradecane standard))	≤ 2	ng/ml
GC/MSD (retention range n-undecane to n-tetracontane, scanning area 30 - 600 amu, individual signals (n- tetradecane standard))	≤ 2	ng/ml

Jeannette David

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.