



Specification

1.16899.0001 Ammonium Test for RQflex® 20
Method: reflectometric with test strips and reagents 5.0 - 20.0 mg/l NH₄⁺
Reflectoquant®

Specification

function control	conforms	
Assay (standard 1)	6.5 - 8.5	mg/l
Assay (standard 2)	9.0 - 13.0	mg/l
Assay (standard 3)	13.5 - 17.5	mg/l

Sharon Brilmayer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.