



# Specification

---

1.19814.0000 Fluoride standard solution traceable to SRM from NIST  
NaF in H<sub>2</sub>O 1000 mg/l F Certipur®

---

## Specification

---

Concentration	$\beta$ (F <sup>-</sup> )	990 - 1010	mg/l
---------------	---------------------------	------------	------

Determination method: Titration with Lanthannitrate solution  
(traceable to NIST SRM 682)  
Accuracy of the method:  $\pm$  2 mg/l

Ayfer Yildirim

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.