



# Specification

---

8.52094.0000 Fmoc-Lys(Mmt)-OH Novabiochem®

---

Specification		
Color (visual)	white to slight yellow to beige	
Appearance of substance (visual)	powder	
Identity (IR)	passes test	
Enantiomeric purity	≥ 99.0	% (a/a)
Purity (TLC(0009))	≥ 95	%
Assay (HPLC, area%)	≥ 95.0	% (a/a)
Assay (acidimetric)	≥ 80.0	%
Water (K. F.)	≤ 2.0	%

For TLC-Systems see our General Catalog.

Product is stabilised with up to 20% diisopropylether upon drying.

Dr. Jörg Bauer

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.