



Specification

8.55005.0000 NovaSyn®TGA resin (90 µm) Novabiochem®

Specification

Color (visual)	white to pale yellow to beige	
Appearance of substance (visual)	beads	
Loading (determined from the substitution of the Fmoc-Leu loaded resin)	0.20 - 0.30	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g

The base resin is PEG-PS-copolymer (90µm) functionalized with 4-hydroxy-methylphenoxyacetic acid.

Dr. Jörg Bauer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.