



# Specification

---

8.55051.0000 Biotin NovaTag™ resin Novabiochem®

---

---

## Specification

---

Color (visual)	white to light yellow to beige	
Appearance of substance (visual)	beads	
Loading (Photometric determination of the Fmoc-chromophore liberated upon treatment with DBU/DMF)	0.35 - 0.65	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g

The polymer matrix is copoly (styrene - 1% DVB), 100 - 200 mesh.

Dr. Jörg Bauer

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.