



# Specification

---

8.55058.0000 Universal PEG NovaTag™ resin Novabiochem®

---

Specification		
Color (visual)	white to light yellow to beige	
Appearance of substance (visual)	beads	
Loading (Photometric determination of the Fmoc-chromophore liberated upon treatment with DBU/DMF)	0.20 - 0.50	mmol/g
Swelling Volume (in DMF)	lot specific result	ml/g

Dr. Jörg Bauer

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.