

Product Specifications

Product name:

Monoclonal Anti-Heat Shock Protein 70 antibody produced in mouse

monoclonal BRM-22, ascites fluid

Product no: **H5147**

Test:

Specifications:

Country of origin

Antibody origin

Israel

Appearance

Form

Liquid

Clarity

Clear

Slightly hazy

Color

Yellow brown

Protein concentration by Biuret

Protein concentration by Biuret

15 - 55 mg/ml

Isotype concentration by RID

Mouse IgG1 concentration

2 - 4 mg/ml

Immunoblotting

Cell tissue extract/Antigen

Bovine brain extract

Second antibody

A9917 Anti Mouse IgG (Fab)-Peroxidase

A2304 Anti Mouse IgG (Fab)-Peroxidase

Substrate

ECL

ECL

Working dilution

1/5000 - 1/10000

MW detected (specific band) approx.

70 kDa

Recommended Retest

Recommended Retest period

5 Years