

## Product Specifications

Product name:

**Monoclonal Anti-Tumor Necrosis Factor- $\alpha$  antibody produced in hamster**

~ 2mg/mL, monoclonal TN3-19.12, purified immunoglobulin, buffered aqueous solution

Product no: **T2824**

### Test:

### Specifications:

#### Country of origin

Antibody origin Israel

#### Appearance

Form Liquid  
Clarity Clear  
Color Colorless  
Odor None

#### Protein Concentration by OD at 280nm

IgG concentration 1.8 - 2.0 mg/ml

#### Immunoblotting (chemiluminescent)

Cell extract/antigen T-7539 recombinant alpha TNF  
recombinant alpha TNF  
Second antibody HRP Mc anti Hamster IgG  
Working concentration 0.25 - 0.5  $\mu$ g/ml  
MW detected (specific band) approx. 17 kDa

#### Recommended Retest

Recommended Retest period 4 Years