


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dichlorodimethylsilane - $\geq 99.5\%$

Product Number: 440272
CAS Number: 75-78-5
MDL: MFCD00000491
Formula: $\text{C}_2\text{H}_6\text{Cl}_2\text{Si}$
Formula Weight: 129.06 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (GC)
Expiration Date Period
5 Years

Colorless to Faint Yellow
Liquid
Conforms to Structure
 $\geq 99.50\%$

Specification: PRD.2.ZQ5.10000012901

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.