

Product Specification

Product Name :

EX-CELL® 293 Serum-Free Medium for HEK 293 Cells Serum-free, Animal-protein free, without L-glutamine, liquid, suitable for cell culture

Product Number: 14571C**MDL Number :** MFCD01634638**Storage Temperature :** +2°C to +8°C

Test	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
Appearance (Form)	Solution
pH	7.0 - 7.4
Osmolality	260 - 300 mOs/kg
Sterility by USP guidelines	No Microbial Growth Detected
Endotoxin Level	≤ 1.0EU/ml
Growth Promotion(Cell Density)	≥ 2000000
Cell Line	HEK 293 Cell Line
Cytotoxicity	Non-toxic

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1266512

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

