

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hybridization Solution – for Northern and Southern blotting, liquid (2x concentrate)

Product Number:

H7140

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

DNase, Exonuclease Detection

None Detected

NICKase, Endonuclease Detection

None Detected

RNAse Detection

None Detected

Suitability

Suitable

for Nucleic Acid Hybridization

Detects at least 50 pg of target DNA

Specification: PRD.2.ZQ5.10000012445

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.