

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Poly(ethylene glycol) dimethacrylate – average M_n 1,000, contains MEHQ as inhibitor

Product Number:

906859



Storage Temperature:

-20 °C

TEST

Specification

Appearance (Color)

White to Off White

Appearance (Form)

Powder or Semi-Solid or Chunk(s)

 ^1H NMR Spectrum

Conforms to Structure

Molecular Number by HNMR

800 - 1200

Expiration Date Period

5 Years

Specification: PRD.0.ZQ5.10000103015

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

