

Product Specification

Product Name : 5R-PLEX beads for PCR purification

Product Number: MBD6009

Storage Temperature : Cooler/Refrigerated

Test	Specification
Suitability	Pass
PCR product clean-up-retaining \geq 200bp DNA fragments	
Appearance (Form)	Suspension
Appearance (Color)	Brown to Very Dark Brown and Very Dark Brown-Red
Recommended Retest Period	-----
2 Years	

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1092224

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

