


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Xylene Cyanol FF - for molecular biology, BioReagent

Product Number: **X4126**
CAS Number: 2650-17-1
MDL: MFCD00019481
Formula: C₂₅H₂₇N₂NaO₆S₂
Formula Weight: 538.61 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Light Green to Very Dark Green and Purple to Very Dark Purple and Black
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Turbidity)	Clear
c = 1mg/ml; Water	
Solubility (Color)	Dark Blue
Suitability	Suitable
Dye Content	
Approximately 75%	
Recommended Retest Period	-----
4 Years	

Specification: PRD.5.ZQ5.10000026801

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

