


Product Specification

Product Name:

Glutaraldehyde solution - technical, ~50% in H₂O (5.6 M)

Product Number:

49629



Formula Weight:

100,12 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid or Viscous Liquid
Appearance (Turbidity)	Clear
Titration	48.0 - 52.0 %
Aldehyde Method	
Refractive index at 20 °C	1.419 - 1.423

Specification: PRE.0.ZQ5.10000024913

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

