

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sigmacote® – siliconizing reagent for glass and other surfaces

Product Number: SL2

Formula Weight: 100.21 g/mol
Storage Temperature: 2 - 8 °C

| TEST | Specification |
|---|---------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Solution |
| Appearance (Turbidity) | Clear |
| Essentially free of particles/precipitate | |
| Performance | Conforms |
| When one drop of an aqueous dye solution is added to glass vial coated with Sigmacote, the dye solution forms a bead and does not stain the vial. | |

Specification: PRD.0.ZQ5.10000034736

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.