

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Peptide bioink - suitable for 3D bioprinting applications

Product Number:

916110

Storage Temperature:

-20 °C

TEST	Specification
Kit Comp (Component 1)	-----
Bioink Powder	
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Kit Comp (Component 2)	-----
Prep Solution	
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Expiration Date Period	-----
16 Months	

Specification: PRD.0.ZQ5.10000104576

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

