

## Product Specification

**Product Name :** Mitomycin C ready made solution 10 mg/mL in DMSO, ≥98% (HPLC), Cell culture tested  
**Product Number:** M5353  
**Storage Temperature :** Freezer -20°C

Test	Specification
Appearance (Form)	Liquid
Appearance (Color)	Light Purple to Very Dark Purple
Purity (HPLC)	≥ 98.0%
Endotoxin Level (LAL)at 1 ml/L ≤0.5 EU/ml	Pass
Handling Precautions Product is light sensitive	Conforms
Cell Culture Test (Suitability)	Pass
Recommended Retest Period 2 Years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.03.ZQ5.1095522

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact [mlsbranding@sial.com](mailto:mlsbranding@sial.com)

