

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

 $\beta$ -Glucanase 2, thermostable – recombinant, expressed in E. coli,  $\geq 90\%$  (SDS-PAGE)

Product Number:

**G8673**

CAS Number:

62213-14-3

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

mg protein/ml (Bradford)

 $\geq 20$ 

un/mg protein (Enzyme Activity)

 $\geq 1.0$ 

One unit will produce 1 umole of reducing sugar (measured as glucose) from Beta-glucan per minute at pH 5.8 at 70 Deg C.

Purity by SDS-PAGE

 $\geq 90\%$ 

Specification: PRD.0.ZQ5.10000049551

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.