

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Liposome Kit: Lipid mixtures for the preparation of liposomes – lyophilized powder

Product Number:

L4395

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Faint Yellow
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
After reconstitution of one vial in 5 mL, CHCl ₃	
Kit Comp (Component 1)	Conforms
L-A-Phosphatidylcholine (egg yolk)	
63 umoles/vial	
Kit Comp (Component 2)	Conforms
Stearylamine	
18 umoles/vial	
Kit Comp (Component 3)	Conforms
Cholesterol	
9 umoles/vial	

Specification: PRD.0.ZQ5.10000090601

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.