

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

SYBR® Green II RNA gel stain – 10,000 × in DMSO

Product Number:

S9305

CAS Number:

172827-25-7

MDL:

MFCD00284716

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

Orange

Appearance (Form)

Solution

Optical Density

0.52 - 0.60

at 481 nm (+ /- 3 nm)

(1206x dilution with 50 mM potassium phosphate buffer, pH 7)

Registered Trademark

SYBR is a registered trademark of Molecular Probes, Inc.

Specification: PRD.0.ZQ5.10000011767

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.