

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trypsin inhibitor from Glycine max (soybean) – lyophilized powder

Product Number:

T9128

CAS Number:

9035-81-8

MDL:

MFCD00082103

Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Color)

White to Light Yellow to Light Brown and
Faint Grey to Light Grey

Appearance (Form)

Powder

mg Trypsin/mg Inhibitor

≥ 1.0

One mg will inhibit a minimum of 1.0 mg Trypsin with an activity of approx. 10,000 BAEE units per mg protein

Recommended Retest Period

6 years

Specification: PRD.2.ZQ5.10000031770

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

