

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium Phosphate Monobasic solution – reagent grade, 1.0 M in solution

Product Number:

P8709 KH_2PO_4

MDL:

MFCD00011401

Formula Weight:

136.09 g/mol

TEST	Specification
------	---------------

Appearance (Color)

Colorless

Appearance (Turbidity)

Clear

Appearance (Form)

Liquid

Concentration

0.97 - 1.03 M

Specification: PRD.0.ZQ5.10000033411

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.