

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

TissueFab® GelAlg – LAP Bioink - low endotoxin, 0.2 µm filtered, suitable for 3D bioprinting applications

Product Number:

925055

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	Colorless to Faint Yellow
Appearance (Form)	Viscous Liquid or Gel
¹ H NMR Spectrum	Conforms to Structure
Bioburden	< 10 CFU/g
Aerobic	
Bioburden	< 10 CFU/g
Fungal	
Endotoxin Level	< 100 EU/ml
Endotoxin units	
Registered Trademark	Note
TissueFab is a registered trademark of Sigma-Aldrich Co. LLC	

Specification: PRD.0.ZQ5.10000107361

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

