

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Ficoll® 400 – lyophilized powder, Type 400-DL, BioReagent, suitable for cell culture

Product Number:

F8016

CAS Number:

26873-85-8

MDL:

MFCD00081599

Formula Weight:

400 kg/mol

### TEST

### Specification

Appearance (Color)

White to Off White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Very Faint Yellow

Solubility (Turbidity)

Clear to Slightly Hazy

100 mg/ml H<sub>2</sub>O

µmhos/cm Conductivity

≤ 200

10% solution

Recommended Retest Period

5 years

Cell Culture Test

Pass

Note

Ficoll is a registered trademark of Cytiva

Specification: PRD.2.ZQ5.10000009668

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

