

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Fluorometric Intracellular Ros Kit – sufficient for 200 fluorometric tests (orange)

Product Number:

MAK144

Storage Temperature:

-20 °C

| TEST | Specification |
|------|---------------|
|------|---------------|

Expiration Date

Confirmed

Report Result

Suitability

Suitable

Suitable for detection of thiols

Specification: PRD.1.ZQ5.10000103499

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

