

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Myelin Basic Protein (MBP) antibody produced in rabbit – ~1 mg/mL, affinity isolated antibody, buffered aqueous solution

Product Number:

M3821

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)

in PBS containing 0.02% sodium azide

Concentration

Western Blot

Immunohistochemistry

Solution

1 - 1 mg/ml

Dilution Range 1:500 - 1:1000

Dilution Range 1:100 - 1:200

Specification: PRD.1.ZQ5.10000009334

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.