

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Gold nanoparticles - 10 nm diameter, OD 1, suspension in 0.1 mM PBS, reactant free

Product Number: 752584

Formula: Au
Formula Weight: 196.97 g/mol
Storage Temperature: 2 - 8 °C

Au

TEST	Specification
Appearance (Form)	Liquid or Suspension
ICP Major Analysis	Confirmed
Confirms Au Component	
Absorption Max	515 - 519 nm
Polydispersity Index (PDI)	≤ 0.2
Core Size	8 - 12 nm
Mean Hydrodynamic Diameter (Z)	11 - 25 nm
Particles/ml	Confirmed
5.38E+ 12 - 6.58E+ 12	
Concentration	Confirmed
(OD): 1	
Buffer	Confirmed
Stabilized suspension in 0.1 mM phosphate buffered saline	
Product of Supplier	Confirmed
CytoDiagnostics, Inc	
Recommended Retest Period	-----
9 months	

Specification: PRD.4.ZQ5.10000100475

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

