

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 γ -Globulins from bovine blood – $\geq 99\%$ (agarose gel electrophoresis)

Product Number:

G5009

CAS Number:

9007-83-4

MDL:

MFCD00163150

Storage Temperature:

-20 °C

TEST**Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Solubility (Color)

Colorless to Faint Yellow

Solubility (Turbidity)

Clear to Hazy

25 mg/ml, 0.9% NaCl

Water (by Karl Fischer)

 $\leq 5\%$

Nitrogen

14 - 16 %

Salt Content

 $\leq 4\%$

Sodium chloride (Based on Na by ICP)

Agarose Electrophoresis

 $\geq 99\%$

Recommended Retest Period

3 years

Specification: PRD.2.ZQ5.10000008824

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.