

Product Specification

Product Name :

Oleic Acid-Albumin from bovine serum liquid, sterile-filtered, BioReagent, suitable for cell culture

Product Number: O3008**MDL Number :** MFCD00284022**Storage Temperature :** +2°C to +8°C

Test	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
pH	6.5 - 7.5
Osmolality Report result	
Endotoxin Level	≤ 400EU/ml
Oleic Acid	≥ 850µg/ml
Protein by Biuret	90 - 120 mg/mL
UV Absorbance 280nm 1:100 dilution	0.60 - 0.79
Sterility by USP guidelines	Negative
Cell Toxicity Assay at 10mL/L	Pass
Chapter 4(D) - Heat Product meets European Union requirements for treated technical blood products. Form of treatment heat treated at 65 deg C for 3 hours.	Conforms

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1267401

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

