

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
pBR322 Plasmid DNA from E. coli RRI - lyophilized powder

Product Number: **D9893**
CAS Number: 93384-17-9
MDL: MFCD00146650

Storage Temperature: 2 - 8 °C

TEST	Specification
A280/A260 UV Absorption	≥ 1.8
Restriction Enzyme Digest	Pass
Complete digestion and correct band pattern observed with EcoR I restriction endonuclease.	

Specification: PRD.1.ZQ5.10000066772

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.