

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N-Hydroxysuccinimidyl-Sepharose® 4 Fast Flow – matrix fast-flow highly cross-linked 4% beaded agarose

Product Number:

H8280

MDL:

MFCD00164950

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Suspension

ml suspension/ml packed gel

1.1 - 1.5

Gel Content

Ligand Immobilized

16 - 23

umol/mL

Note

Sepharose is a registered trademark of Cytiva

Specification: PRD.3.ZQ5.10000008658

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

