

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

TRI Reagent® - For processing tissues, cells cultured in monolayer or cell pellets

Product Number:

T9424

MDL:

MFCD00213058

TEST	Specification
------	---------------

Appearance (Color)

Red

Appearance (Form)

Liquid

Ratio

≥ 1.70

Spectrophotometric Analysis of RNA

260/280 Ratio

Bioassay

Conforms

Bioanalyzer Analysis: Shows two intact ribosomal bands.

Registered Trademark

Confirmed

TRI Reagent is a Registered Trademark of Molecular Research Center, Inc.

Specification: PRD.3.ZQ5.10000004745

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.