

Product Specification

Product Name : Anti-CD4 antibody, Mouse monoclonal clone Q4120, purified from hybridoma cell culture
Product Number: C1805
Storage Temperature : Cooler/Refrigerated

Test	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	50 - 200 µg/ml
FACS -Cells	Human periferal blood lymphocytes
FACS-Working Volume µl/tst	5 - 5
Recommended Retest Period 3 Years	-----

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.04.ZQ5.1106923

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

