

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

NanoFabTX™ -DOTAP Lipid Mix - for synthesis of cationic (DOTAP) liposomes

Product Number:

926027

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Paste or Solid
Recommended Retest Period	-----
15 months	
Trademark	Note
NanoFabTx is a trademark of Sigma-Aldrich Co. LLC	

Specification: PRD.0.ZQ5.10000108102

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

