

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Glutamate Receptor NMDAR2C (NR2C) antibody produced in rabbit – affinity isolated antibody, lyophilized powder

Product Number: M266

MDL: MFCD01865098

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Lyophilized from 5 mM ammonium bicarbonate.

Working Dilution	Conforms
------------------	----------

Immunoblotting, Immunohistochemistry and immunoprecipitation: 1:1,000

Immunoprecipitation	Conforms
---------------------	----------

Recommended working dilution: 3 ul per 200 ug lysate.

Specification: PRD.0.ZQ5.10000055783

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.