

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
HIS-Select® Nickel Affinity Gel

Product Number: **P6611**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	Faint Blue to Very Dark Blue
Appearance (Form)	Suspension
μmol Ni/mL packed gel	6 - 18
Binding Capacity / ml PGV	≥ 15
> or = 15 mg target protein/mL packed gel	
ml suspension/ml packed gel	1.5 - 2.4
(mL resin slurry will give 1 mL pgv),	
Settling Volume (Resin Slurry)	

Specification: PRD.2.ZQ5.10000006060

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.