


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SYLGARD® 184 - 10 g clip-pack

Product Number: 761036



TEST	Specification
Appearance (Curing Agent) Colorless to Light Straw	Pass
Appearance (Base) Colorless to Light Straw Liquid	Pass
Viscosity (cst)	4000.00 - 6500.00
Registered Trademark Sylgard(R) 184 Silicone Elastomer Kit is a Registered Trademark of Dow Corning Corporation	Confirmed
Note All test results are in regards to the Base unless otherwise stated.	Confirmed
Expiration Date Period 2 Years from Date of Manufacture	-----

Specification: PRD.2.ZQ5.10000030541

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.